

NEIL SQUIRE SOCIETY'S



GameStik Setup Guide

Version 1.11 (October 20,2018)



Table of Contents NEIL SQUIRE SOCIETY'S



 	1			
SameStik Setup Guide				
Required Software:				
Installation Instructions:				
Functions (Computer):				
Functions (Xbox Adaptive Controller):				
References	5			



Required Software:

GameStik firmware

https://github.com/makersmakingchange/GameStik/blob/master/Software/GameStik_Firmware/GameStik Firmware.ino

The Joystick library needs to be installed before uploading Joystick firmware from following link:

https://github.com/MHeironimus/ArduinoJoystickLibrary

Installation Instructions:

1. Download the GameStik firmware from following link:

https://github.com/makersmakingchange/GameStik/blob/master/Software/GameStik_Firmware/GameStik Firmware.ino

2. Visit the "ArduinoJoystickLibrary" GitHub repository using following link: [1]

https://github.com/MHeironimus/ArduinoJoystickLibrary

- Download the "ArduinoJoystickLibrary" by clicking on "Clone or Download" button and then selecting "Download ZIP" option.
- 4. Extract the "ArduinoJoystickLibrary-master.zip" file.
- 5. Copy the Joystick folder to the Arduino libraries folder.

An example of path to Arduino libraries folder in Windows: "C:\Program Files (x86)\Arduino\libraries".

- 6. Once the folder has been copied, the Joystick library should appear in the Arduino IDE list of libraries.
- 7. Open Arduino IDE and you should be able to see Joystick examples under "File>Examples>Joystick".
- 8. Now open the "GameStik_Firmware.ino" file using Arduino IDE.
- Select the "Arduino/Genuino Micro" as the Board and then select the available port number under Tools option.
- 10. Click on upload button to upload the code.



Functions (Computer):

Action	Function	Time (Seconds)
Joystick (Right)	Move Right on X-axis	N/A
Joystick (Down)	Move Down on Y-axis	N/A
Joystick (Left)	Move Left on X-axis	N/A
Joystick (Up)	Move Up on Y-axis	N/A
Puff	Press Button Number 5	1
Sip	Press Button Number 6	1
Puff	Press Button Number 4	2
Sip	Press Button Number 1 ¹	2
Puff	Press Button Number 7	3
Sip	Press Button Number 8	3

Table 1: Computer Functions

Note 1: The press button 1 or 2-second long sip is reserved to be used for Xbox Adaptive Controller.



Functions (Xbox Adaptive Controller):

The following functions are the default options for the specified actions which can be modified using profile settings. [2]

Action	Function (Left USB)	Function (Right USB)	Time (Seconds)
Joystick (Right)	Move Right on X-axis	Move Right on X-axis	N/A
Joystick (Down)	Move Down on Y-axis	Move Down on Y-axis	N/A
Joystick (Left)	Move Left on X-axis	Move Left on X-axis	N/A
Joystick (Up)	Move Up on Y-axis	Move Up on Y-axis	N/A
Puff	А	X	1
Sip	В	Y	1
Puff	LB	RB	2
Sip	X1 ²	View ¹	2
Puff	View	X1	3
Sip	Menu	X2	3

Table 2: Xbox Adaptive Controller Functions

The Xbox Adaptive Controller comes with the option to configure and map buttons using three different mapping profiles. The profiles can be switched using the *Profile* button on the Xbox Adaptive Controller. The profiles can be edited based on the needs of the user. [3] The *Shift* button action can be used to switch between buttons when the *Shift* functionality is configured in the *Profile* settings. [4]

Note 2: The LED on the GameStik will blink once in green when the X1 button action or 2-second long sip is performed. The X1 button action or 2-second long sip is reserved to be used as Shift button which is used to switch between primary and shift button actions.

Please follow the instructions below to map the buttons through the *Profile* settings:

- Press the Xbox button to open the guide.
- Select Home > Settings > Kinect and Devices and choose Devices & Accessories.
- Press the A button or short puff.
- Select Xbox Adaptive Controller.
- Select *Configure* option and then press the *A* button or short puff to enter.
- Select New profile and then press the A button or short puff to enter.



- Enter the name for the new profile. *GameStik1* can be used as an example.
- Select the button that you would like to map under *Mapping* menu and then press the *A* button or short puff to enter.
- Assign the *Primary* and *Shift* button actions for selected button.
- Press the B button or short sip when you are finished.
- The following three examples can be used to map the main buttons:
 - A button or short puff
 - Select the A button under Mapping menu and then press the A button or short puff
 to enter. We select the A button which is activated by short puff as the first example
 when GameStik is connected to the left USB 2.0 port.
 - 2. Assign the *Primary* and the *Shift* button actions for selected *A* button.
 - **3.** The *Primary* button action is assigned to the *A* button and the *Shift* button action is assigned to the *X* button for the *A* button in the *GameStik1* profile example.
 - **4.** Press the **B** button or short sip when you are finished.
 - **B** button or short sip
 - Select the B button under Mapping menu and then press the A button or short puff
 to enter. We select the B button which is activated by short sip as the second
 example when GameStik is connected to the left USB 2.0 port.
 - 2. Assign the *Primary* and the *Shift* button actions for selected *B* button.
 - **3.** The *Primary* button action is assigned to the *B* button and the *Shift* button action is assigned to the *Y* button for the *B* button in the *GameStik1* profile example.
 - **4.** Press the **B** button or short sip when you are finished.
 - LB button or 3 seconds long puff
 - Select the LB button under Mapping menu and then press the A button or short
 puff to enter. We select the LB button which is activated by 2-second long puff as
 the third example when GameStik is connected to the left USB 2.0 port.
 - 1. Assign the *Primary* and the *Shift* button actions for selected *LB* button.
 - The *Primary* button action is assigned to the *LB* button and the *Shift* button action is assigned to the *RB* button for the *LB* button in the *GameStik1* profile example.
 - **3.** Press the **B** button or short sip when you are finished.



- Select the button that you would like to map to the Shift button action under the Mapping menu
 and then press the A button or short puff to enter. We select the X1 button which is activated by 2second long sip as an example when GameStik is connected to the left USB 2.0 port.
- Select the *Use as Shift* button check box and then press the *A* button or short puff to enter.
- Press the **B** button or short sip when you are finished.
- Press the **B** button or short sip again to enter the **Profiles** page again.
- Select the new created profile which is *GameStik1* as an example and then press the *A* button or short puff to enter.
- Next save the new profile by selecting the *No Slot* dropdown menu and then pressing the *A* button
 or short puff to enter.
- Select one of the three available slots and then press the A button or short puff to save the profile.
- GameStik1 profile is now ready to use. Press the Profile button on the Xbox Adaptive Controller to switch between profiles and select GameStik1 profile.

The Table below represents the mapped functions for *GameStik1* profile example.

Action	Primary	Shift	Time (Seconds)
Joystick (Right)	Move Right on X-axis	Move Right on X-axis	N/A
Joystick (Down)	Move Down on Y-axis	Move Down on Y-axis	N/A
Joystick (Left)	Move Left on X-axis	Move Left on X-axis	N/A
Joystick (Up)	Move Up on Y-axis	Move Up on Y-axis	N/A
Puff	А	X	1
Sip	В	Y	1
Puff	LB	RB	2
Sip	X1 (Shift Button)	X1 (Shift Button)	2
Puff	View	View	3
Sip	Menu	Menu	3

Table 3: GameStik1 profile Functions



References

- 1. ArduinoJoystickLibrary Github repository
- 2. What external devices work with the Xbox Adaptive Controller?
- 3. Customize the Xbox Adaptive Controller in the Xbox Accessories app
- **4.** How do I use Shift with the Xbox Adaptive Controller?